

Human HS3ST1 protein (Recombinant) (His-Tag) (STJP018423)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 40% Glycerol, 2mM DTT
Molecular Weight	36.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9957
Gene Symbol	HS3ST1
UniProt ID	HS3S1_HUMAN
Immunogen Region	21-307aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSRPAELGQ QELLRKAGTL QDDVDRDGVAP NGSAAQLPQT IIIGVRKGGT RALLEMLSLH PDVAAAENEV HFFDWEEHYS HGLGWYLSQM PFSWPHQLTV EKTPAYFTSP KVPERVYSMN PSIRLLILR DPSEVLSDY TQVFYNHMQK HKPYPSIEEF LVRDGRNLVD YKALNRSLYH VHMQNWLRFPLRHHIIVDGDRLIRDPF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human HS3ST1 protein 3 µg Coomassie blue stain