

## Human EBNA1BP2 protein (Recombinant) (His-Tag) (STJP018432)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Proteins           |
| Applications  | SDS-PAGE/Denatured |
| Host / Source | E.coli             |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Concentration       | 0.5 mg/mL                                                                                                                                           |
| Formulation         | Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.4M urea                                                                          |
| Molecular Weight    | 37.2 kDa                                                                                                                                            |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                                                |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">10969</a>                                                                                                                                                                                                                                          |
| Gene Symbol        | <a href="#">EBNA1BP2</a>                                                                                                                                                                                                                                       |
| UniProt ID         | <a href="#">EBP2_HUMAN</a>                                                                                                                                                                                                                                     |
| Immunogen Region   | 1-306aa                                                                                                                                                                                                                                                        |
| Immunogen Sequence | MGSSHHHHHH SSGLVPRGSH MGSMDTPPLS DSESESDSL VTDRELQDAF SRGLLKPGLN<br>VVLEGPKKAV NDVNLKQCL AEFKRDLEWV ERLDVTLPV PEIGGSEAPA PQNKDQKAVD PEDDFQREMS<br>FYRQAQAAVL AVLPRLHQLK VPTKRPTDYF AEMAKSDLQM QKIRQLQTK QAAMERSEKA KQLRALRKYG<br>KKVQTEVLQK RQKEAHMMN AIKKYQKG |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



SDS-PAGE gel of recombinant human EBNA1BP2 protein 3ug under reducing conditions