

Human RPL35A protein (Recombinant) (His-Tag) (STJP018438)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	14.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6165
Gene Symbol	RPL35A
UniProt ID	RL35A_HUMAN
Immunogen Region	1-110aa
Immunogen Sequence	MSGRLWS KAIFAGYKRG LRNQREHTAL LKIEGVYARD ETEFYLGKRC AYVYKAKNNT VTPGGKPNKT RVIWGKVTRA HGNSGMVRAK FRSNLPKAI GHRIRVMLYP SRI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant human RPL35A protein SDS-PAGE gel under reducing conditions stained with Coomassie blue