

Human RPL34 protein (Recombinant) (His-Tag) (STJP018439)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M NaCl, 10% Glycerol
Molecular Weight	15.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6164
Gene Symbol	RPL34
UniProt ID	RL34_HUMAN
Immunogen Region	1-117aa
Immunogen Sequence	MVQRLTY RRRLSYNTAS NKTRLSRTPG NRIVLYTKK VGKAPKSACG VCPGRLRGVR AVRPKVLMRL SKTKKHVSRA YGGSMAKCV RDRIKRAFLI EEQKIVVKVL KAQAQSQKAK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human RPL34 protein 3µg under reducing conditions Coomassie stain