

Human DBR1 protein (Recombinant) (His-Tag) (STJP018443)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	38 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51163
Gene Symbol	DBR1
UniProt ID	DBR1_HUMAN
Immunogen Region	228-544aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGS HVKFAAL MQHQA KDKGQ TARATKFLAL DKCLPHRD FL QILEIEHDPS APDYLEYDIE WLTILRATDD LINVTGRLWN MPENNGLHAR WDYSATEEGM KEVLEKLNHD LKVP CNFSVT AACYDPSKPQ TQMQLIHRIN PQTTEFCAQL GIIDINVRLQ KSKEEHVCG EYEEQDDVES NDSGEDQSEY NTDTSALSSI NPDEIMLD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human DBR1 protein under reducing conditions stained with Coomassie blue