

Human Syntaxin 6 protein (Recombinant) (His-Tag) (STJP018446)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% Glycerol, 2mM DTT
Molecular Weight	29.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10228
Gene Symbol	STX6
UniProt ID	STX6_HUMAN
Immunogen Region	1-234aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MGSMSMEDPF FVVKGEVQKA VNTAQGLFQR WTELLQDPST ATREEIDWTT NELRNNLRSI EWDLEDLDET ISIVEANPRK FNLDATELSI RKAFITSTRQ VVRDMKQDMS TSSVQALAEK KNRQALLGDS GSNWSTGTT DKYGRDLREL QRANSHFIEE QQAQQQLIVE QQDEQLELVS GSIGVLKNMS QRIGGELEEQ AVMLEDFS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human Syntaxin 6 protein under reducing conditions stained with Coomassie blue