

Human eIF-3 epsilon/EIF3F protein (Recombinant) (His-Tag) (STJP018460)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.4M urea
Molecular Weight	40kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8665
Gene Symbol	EIF3F
UniProt ID	EIF3F_HUMAN
Immunogen Region	1-357aa
Immunogen Sequence	MATPAVP VSAPPATPTP VPAAAPASVP APTPAPAAAP VPAAAPASSS DPAAAAAATA APGQTPASAQ APAQTPAPAL PGPALPGPFPG GRVVRLHPV ILASIVDSYE RRNEGAARVI GTLLGTVDKH SVEVTNCFVS PHNESEDEVA VDMEFAKNMY ELHKKVSPNE LILGWYATGH DITEHSVLH EYYSREAPNP IHLTVD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human eIF-3 epsilon EIF3F protein STJP018460 under reducing conditions