

Human UBC3B/UBE2R2 protein (Recombinant) (His-Tag) (STJP018477)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	29.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54926
Gene Symbol	UBE2R2
UniProt ID	UB2R2_HUMAN
Immunogen Region	1-238aa
Immunogen Sequence	MAQQQMT SSQKALMLEL KSLQEEPVEG FRITLVDESD LYNWEVAIFG PPNTLYEGGY FKAHIKFPID YPYSPPTFRF LTKMWHPNYI ENGDVCISIL HPPVDDPQSG ELPSEWRNPT QNVRTILLSV ISLLNEPNTF SPANVDASVM FRKWRDSKGG DKEYAEIIRK QVSATKAEAE KDGVKVPTTL AEYCIKTKVP SNDNSS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human UBC3B UBE2R2 protein under reducing conditions Coomassie stain