

Human CHCHD3 protein (Recombinant) (His-Tag) (STJP018480)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% Glycerol, 2mM DTT
Molecular Weight	28.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54927
Gene Symbol	CHCHD3
UniProt ID	MIC19_HUMAN
Immunogen Region	1-227aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMGGTTST RRVTFEADEN ENITVVKGIR LSENVIDRMK ESSPSGSKSQ RYSGAYGASV SDEELKRRVA EELALEQAKK ESEDQKRLKQ AKELDRERAA ANEQLTRAIL RERICSEEEER AKAKHLARQL EEKDRVLKKQ DAFYKEQLAR LEERSSEFYR VTTEQYQKAA EEVEAKFKRY ESHPCADLQ AKILQCYREN THQTLKCS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human CHCHD3 protein under reducing conditions
Coomassie blue stain