

Human RPS18 protein (Recombinant) (His-Tag) (STJP018482)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	20.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6222
Gene Symbol	RPS18
UniProt ID	RS18_HUMAN
Immunogen Region	1-152aa
Immunogen Sequence	MSLVIPE KFAQHILRVLN TNIDGRRKIA FAITAIGVVG RRYAHVVLK ADIDLTKRAG ELTEDEVERV ITIMQNPRQY KIPDWFNLRQ KDVKDGKYSQ VLANGLDNKL REDLERLKKI RAHRGLRHFV GLRVRGQHTK TTGRRGRVTVG VSKKK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RPS18 protein 3µg under reducing conditions Coomassie blue stain