

Human APOBEC4 protein (Recombinant) (His-Tag) (STJP018493)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 0.4M urea
Molecular Weight	44 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	403314
Gene Symbol	APOBEC4
UniProt ID	ABEC4_HUMAN
Immunogen Region	1-367aa
Immunogen Sequence	MEPIYEE YLANHGTIVK PYYWLSFSLD CSNCPYHIRT GEEARVSLTE FCQIFGFYPG TTFPQTKHLT FYELKTSSGS LVQKGHASSC TGNYPHPESM LFEMNGYLDS AIYNNDISIRH IILYSNNSPC NEANHCCISK MYNFLTYPG ITLSIYFSQL YHTEMDFPAS AWNREALRSL ASLWPRVVLV PISGGIWHSV LHSFIS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human APOBEC4 protein under reducing conditions stained with Coomassie blue