

Human eIF-3 beta/EIF3I protein (Recombinant) (His-Tag) (STJP018494)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 0.4M urea
Molecular Weight	38.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8668
Gene Symbol	EIF3I
UniProt ID	EIF3I_HUMAN
Immunogen Region	1-325aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MSGMKPILLQ GHERSITQIK YNREGDLLFT VAKDPIVNVW YSVNGERLGT YMGHTGAVWC VDADWDTKHV LTGSADNSCR LWDCEGTGKQL ALLKTNSAVR TCGDFGGNI IMFSTDKQMG YQCFVSFFDL RDPSQIDNNE PYMKIPCNDK KITSVAVWGPI GECIIAGHES GELNQYSAKS GEVLNVNKEH SRQINDIQLS RDMTMFVT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human eIF-3 beta/EIF3I protein under reducing conditions