

Human MAP1/MOAP1 protein (Recombinant) (His-Tag) (STJP018505)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 0.4M urea
Molecular Weight	41.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64112
Gene Symbol	MOAP1
UniProt ID	MOAP1_HUMAN
Immunogen Region	1-351aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MGSMTLRLLE DWCRGMDMNP RKALLIAGIS QSCSVAEIEE ALQAGLAPLG EYRLGMRMR RDENRKVALV GLTAETSHAL VPKEIPGKGG IWRVIFKPPD PDNTFLSRLN EFLAGEGMTV GELSRALGHE NGSLDPEQGM IPEMWAPMLA QALEALQPAL QCLKYKKLRV FSGRESPEPG EEEFGRWMFH TTQMIKAWQV PDVEKRRR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human MAP1
MOAP1 protein under reducing conditions