

Human MAGEB10 protein (Recombinant) (His-Tag) (STJP018509)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT, 0.1mM PMSF
Molecular Weight	41.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	139422
Gene Symbol	MAGEB10
UniProt ID	MAGBA_HUMAN
Immunogen Region	1-347aa
Immunogen Sequence	MPRGQKS KLRAREKRRQ ARGGLEDLID ALDILEEEEE SPPSASACKL DVFQSSLDGA SNNPHGLREA QSTSTSATAA SHTRHPEGVN DQMEERPCT QDLEATDSFP RGPVDEKVII LVHYLLYKYQ MKEPITKADM LRNVTQMSKS QFPVILSRAS EHLELIFGLD LKEVEPNKHI YVLVNKLDLG CDAKLSDETG VPKTGL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human MAGEB10 protein stained with Coomassie blue