

Human Snail protein (Recombinant) (His-Tag) (STJP018513)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 0.4M urea
Molecular Weight	31.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6615
Gene Symbol	SNAI1
UniProt ID	SNAI1_HUMAN
Immunogen Region	1-264aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMPRSLV RKPSDPNRKP NYSELQDSNP EFTFQQPYDQ AHLAAIPPP EILNPTASLP MLIWDSVLAP QAQPIAWASL RLQESPRVAE LTSLSDEDSDG KGSQPPSPPS PAPSSFSSSTS VSSLEAEAYA AFPGLGQVPK QLAQLSEAKD LQARKAFNCK YCNKEYLSLG ALKMHIRSH LPCVCGTCGK AFSRPWLLQG HVRTHTGE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human Snail protein under reducing conditions Coomassie stain