

Human DcR3/TNFRSF6B protein (Recombinant) (His-Tag) (STJP018516)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 0.4M urea
Molecular Weight	32.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8771
Gene Symbol	TNFRSF6B
UniProt ID	TNF6B_HUMAN
Immunogen Region	30-300aa
Immunogen Sequence	VAETPTY PWRDAETGER LVCAQCPPGT FVQRPCRRDS PTCGPCPPR HYTQFWNYLE RCRYCNVLCG EREEEARACH ATHNRACRCR TGFFAHAGFC LEHASCPPGA GVIAPGTPSQ NTQCQPCPPG TFSASSSSSE QCQPHRNCTA LGLALNVP GS SSHDTLCTSC TGFPLSTRVP GAEECERAVI DFVAFQDISI KRLQRL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing 3µg recombinant human DcR3 TNFRSF6B protein under reducing conditions Coomassie blue stain