

Human RNF34 protein (Recombinant) (His-Tag) (STJP018519)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 0.1M NaCl, 1mM DTT
Molecular Weight	44.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80196
Gene Symbol	RNF34
UniProt ID	RNF34_HUMAN
Immunogen Region	1-373aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMRKAGAT SMWASCCGLL NEVMGTGAVR GQQSAGAGAT GPFRTNPNE FSTYPPAATE GPNIVCKACG LSFVFRKKH VCCDCKKDFC SVCSVLQENL RRCSTCHLLQ ETAFQRPQLM RLKVKDLRQY LILRNIPIDT CREKEDLVDL VLCHHGLGSE DDMDTSSLNS SRSQTSSFFT RSFFSNYTAP SATMSSFQGE LMDGDQTS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human RNF34 protein under reducing conditions Coomassie blue stain