

Human RPL31 protein (Recombinant) (His-Tag) (STJP018533)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% Glycerol, 1mM DTT
Molecular Weight	16.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6160
Gene Symbol	RPL31
UniProt ID	RL31_HUMAN
Immunogen Region	1-125aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMAPAKKG GEKKKGRSAI NEVVTREYTI NIKKRIHGVG FKKRAPRALK EIRKFAMKEM GTPDVIRIDTR LNKAVWAKGI RNVPYRIRVR LSRKRNEDED SPNKLYTLVT YVPVTTFKNL QTVNVN DEN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RPL31 protein under reducing conditions stained with Coomassie blue