

Human BLC/BCA-1/CXCL13 protein (Recombinant) (His-Tag) (STJP018537)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing, 10% Glycerol, 0.4M urea
Molecular Weight	12.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10563
Gene Symbol	CXCL13
UniProt ID	CXL13_HUMAN
Immunogen Region	23-109aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSVLEVYYT SLRCRCVQES SVFIPRRFID RIQILPRGNG CPRKEIIVWK KNKSIVCVDV QAEWQRMME VLRKRSSSTL PVPVFKRKIP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human BLC BCA-1 CXCL13 protein under reducing conditions