

Human BEND6 protein (Recombinant) (His-Tag) (STJP018541)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.1M NaCl
Molecular Weight	33.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	221336
Gene Symbol	BEND6
UniProt ID	BEND6_HUMAN
Immunogen Region	1-279aa
Immunogen Sequence	MQKIVQT DEITNTQAFR KGKRKRTETM DSENANSMDMD KGQRDPYSGN AFLPGESSSE DEEPLAELSK EELCAKIKSL KEKLTNTRKE NSRLRQSLVM LQVLPQAVTQ FEELVGMAEA LLKGGGTMTS SASTLWRATN NSSPDSFAST CSNSNSNSSS PVSLKPEEEH QTDEKQFQIE KWQIARCNKS KPQKFINLDM QVLYTN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human BEND6 protein 3 µg under reducing conditions
Coomassie blue stain