

Human FAM84B/LRATD2 protein (Recombinant) (His-Tag) (STJP018542)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.15M NaCl, 1mM DTT, 0.1mM PMSF
Molecular Weight	36.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	157638
Gene Symbol	LRATD2
UniProt ID	LRAT2_HUMAN
Immunogen Region	1-310aa
Immunogen Sequence	MGNQVEK LTHLSYKEVP TADPTGVDRD DGPRIGVSYI FSNDDDEDVEP QPPPQGPDDG GLPDGGDGGP PPQPQPYDPR LHEVECSVFY RDECIYQKSF APGSAALSTY TPENLLNKCK PGDLVEFVSQ AQYPHWAVVYV GNFQVVHLHR LEVINSFLTD ASQGRRRGRVV NDLYRYKPLS SSAVVRNALA HVGAKERELS WRNSES

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human FAM84B LRATD2 protein 3µg under reducing conditions with Coomassie blue stain