

Human Dynactin Subunit 2/DCTN2 protein (Recombinant) (His-Tag) (STJP018550)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	46.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10540
Gene Symbol	DCTN2
UniProt ID	DCTN2_HUMAN
Immunogen Region	1-403aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MGSMADPKYA DLPGIARNEP DYYETSDLPE DDQAEFDAEL EELTSTSVEH IIVNPNAAYD KFKDKRVGK GLDFSDRIGK TKRTGYESGE YEMLGEGLGV KETPQQKYQR LLHEVQELTT EVEKIKTTVK ESATEEKLTP VLLAKQLAAL KQLVASHLE KLLGPDAAIN LTDPDGALAK RLLQLEATK NSKGGSGGKT TGTPPDSS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human Dynactin Subunit 2 DCTN2 protein under reducing conditions