

Human NRAGE protein (Recombinant) (His-Tag) (STJP018555)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 1mM DTT
Molecular Weight	31.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9500
Gene Symbol	MAGED1
UniProt ID	MAGD1_HUMAN
Immunogen Region	504-760aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSLRPSPNS RASQNPAAQ PRDVALLQER ANKLVKYLML KDYTKVPIKR SEMLRDIIIE YTDVYPEIIE RACFVLEKKF GIQLKEIDKE EHLIYLISTP ESLAGILGTT KDTPKLGLLL VILGVIFMNG NRASEAVLWE ALRKMGLRPG VRHPLLGDLR KLLTYEFVKQ KYLDYRRVNP SNPPEYEFW GLRSYHETSK MKVLRFA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human NRAGE protein under reducing conditions with Coomassie blue stain