

Human RAD51B protein (Recombinant) (His-Tag) (STJP018565)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing, 10% Glycerol, 0.4M urea
Molecular Weight	40.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5890
Gene Symbol	RAD51B
UniProt ID	RA51B_HUMAN
Immunogen Region	1-350aa
Immunogen Sequence	MGSKKLK RVGLSQELCD RLSRHQILTC QDFLCCLSPLE LMKVTGLSYR GVHELLCMVS RACAPKMQTA YGIKAQRSAD FSPAFLSTTL SALDEALHGG VACGSLTEIT GPPGCGKTQF CIMMSILATL PTNMGGLEGA VVYIDTESAF SAERLVEIAE SRFPRYFNTE EKLLLTSSKV HLYRELTCD E VLQRIESLEE EIISKG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human RAD51B protein 3ug under reducing conditions with Coomassie blue stain