

Human Apolipoprotein A5/APOA5 protein (Recombinant) (His-Tag) (STJP018576)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing, 10% Glycerol, 0.4M urea
Molecular Weight	41.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	116519
Gene Symbol	APOA5
UniProt ID	APOA5_HUMAN
Immunogen Region	24-366aa
Immunogen Sequence	RKGFWDY FSQTSGDKGR VEQIHQQKMA REPATLKDSL EQDLNNMNKF LEKLRPLSGS EAPRLPQDPV GMRRQLQEEL EEVKARLQPY MAEAHELVGW NLEGLRQQLK PYTMDLMEQV ALRVQELQEQ LRVVGEDTKA QLLGGVDEAW ALLQGLQSRV VHHTGRFKEL FHPYAESLVS GIGRHVQELH RSVAPHAPAS PARLSR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human APOA5 protein under reducing conditions stained with Coomassie blue