

Human JAK2 protein (Recombinant) (His-Tag) (STJP018587)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing, 10% Glycerol, 0.4M urea
Molecular Weight	18.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3717
Gene Symbol	JAK2
UniProt ID	JAK2_HUMAN
Immunogen Region	1014-1132aa
Immunogen Sequence	MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSHMGE SPIFWYAPES LTESKFVSAS DVWSFGVVLY ELFTYIEKSK SPPAEFMRMI GNDKQGQMIV FHLLIELLKNN GRLPRPDGCP DEIYMIMTEC WNNNVNQRPS FRDLALRVDQ IRDNMAG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human JAK2 protein under reducing conditions Coomassie blue stain