

Human Dynactin Subunit 2/DCTN2 protein (Recombinant) (His-Tag) (STJP018590)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 0.15M NaCl, 1mM DTT
Molecular Weight	47.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10540
Gene Symbol	DCTN2
UniProt ID	DCTN2_HUMAN
Immunogen Region	1-406aa
Immunogen Sequence	MADPKYA DLPGIARNEP DVIYETSDLPE DDQAEFDFA QELEELTSTS VEHIIVNPNA AYDKFKDKRV GTKGLDFSDR IGKTKRTGYE SGEYEMLGEG LGVKETPQQK YQRLLHEVQE LTTEVEKIKT TVKESATEEK LTPVLLAKQL AALKQQLVAS HLEKLLGPDA AINLTDPDGA LAKRLLQLQ ATKNSKGGSG GKTTGT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Dynactin Subunit 2 DCTN2 protein 3µg Coomassie stain