

Human Galectin-13/LGALS13 protein (Recombinant) (His-Tag) (STJP018593)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 30% Glycerol, 1mM DTT
Molecular Weight	18.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29124
Gene Symbol	LGALS13
UniProt ID	PP13_HUMAN
Immunogen Region	1-139aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMSLPVP YKLPVLSVG SCVHIKGTPI HSFINDPQLQ VDFYTDMDDED SDIAFRFRVH FGNHVVMMNRR EFGIWMLEET TDYVPFEDGK QFELCIYVHY NEYEIKVNGI RIYGFVHRIP PSFVKMVQVS RDISLTSVCV CN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing recombinant human Galectin-13 LGALS13 protein 3µg under reducing conditions