

Human CEND1 protein (Recombinant) (His-Tag) (STJP018594)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	15kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51286
Gene Symbol	CEND1
UniProt ID	CEND_HUMAN
Immunogen Region	1-125aa
Immunogen Sequence	MESRGKS ASSPKPDTKV PQVTTEAKVP PAADGKAPLT KPSKKEAPE KQPPAAPTT APAKKTSAKA DPALLNNHSN LKPAPTVPS PDATPEPKGP GDGAEDEEA SGGPGGRGPW SCENFNPL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human CEND1 protein 3µg under reducing conditions with Coomassie blue stain