

Human MLF1IP protein (Recombinant) (His-Tag) (STJP018595)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.4M urea
Molecular Weight	33.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79682
Gene Symbol	CENPU
UniProt ID	CENPU_HUMAN
Immunogen Region	147-418a
Immunogen Sequence	TRRKVKS AEKISTQRHE VIRTASSEK SEKPAESVTS KKTGPLSAQP SVEKENLAIE SQSKTQKKGK ISHDKRKKSR SKAIGSDTSD IVHIWCPEGM KTSIDIKELNI VLPEFEKTHL EHQQRIESKV CKAAIATFYV NVKEQFIKML KESQMLTNLK RKNAMISDI EKKRQRMIEV QDELLRLEPQ LKQLQTKYDE LKERKS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human MLF1IP protein under reducing conditions stained with Coomassie blue