

Human BRMS1 protein (Recombinant) (T7-Tag) (STJP018602)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.4M urea
Molecular Weight	29.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25855
Gene Symbol	BRMS1
UniProt ID	BRMS1_HUMAN
Immunogen Region	1-246aa
Immunogen Sequence	MASMTGGQQM GRGSHMPVQP PSKDTEEMEA EGDSAAEMNG EEEEEEEERS GSQTESEES SEMDDDEDYER RRSECVSEML DLEKQFSELK EKLFRERLSQ LRLRLEEVG ERAPEYTEPL GGLQRSLKIR IQVAGIYKGF CLDVIRNKYE CELQGAQHL ESEKLLYDT LQELQERIQ RLEEDRQSLD LSSEWWDDKL HARGSSRSWD SLPPSKRKKA PLVSGPYI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human BRMS1 protein under reducing conditions Coomassie blue stain