

Human SORBS3 protein (Recombinant) (His-Tag) (STJP018604)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 1mM DTT
Molecular Weight	39.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10174
Gene Symbol	SORBS3
UniProt ID	VINEX_HUMAN
Immunogen Region	1-329aa
Immunogen Sequence	MADGGSP FLGRRDFVYP SSTRDPSASN GGGSPARREE KKRKAARLKF DFQAQSPKEL TLQKGDIVYI HKEVDKNWLE GEHHGRLGIF PANYVEVLP DEIPKPIKPP TYQVLEYGEA VAQYTFKGD L EVELSFRKGE HICLIRKVNE NWYEGRITGT GRQGI FPASY VQVSREPRLR LCDDGPQLPT SPRLTAAARS ARHPSS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human SORBS3 protein stained with Coomassie blue under reducing conditions