

Human GNG4 protein (Recombinant) (His-Tag) (STJP018619)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	10.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2786
Gene Symbol	GNG4
UniProt ID	GBG4_HUMAN
Immunogen Region	1-72aa
Immunogen Sequence	MGSSHHHHH SSSLVPRGSH MGSMKEGMSN NSTTISQAR KAVEQLKMEA CMDRVKVSQA AADLLAYCEA HVREDPLIIP VPAENPFRE KKFFC

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human GNG4 protein under reducing conditions stained with Coomassie blue