

## Human MMADHC protein (Recombinant) (His-Tag) (STJP018629)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | SDS-PAGE |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Concentration       | 1 mg/mL                                                                                                                                             |
| Formulation         | Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT                                                                |
| Molecular Weight    | 31 kDa                                                                                                                                              |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                                                         |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">27249</a>                                                                                                                                                                                                                                                   |
| Gene Symbol        | <a href="#">MMADHC</a>                                                                                                                                                                                                                                                  |
| UniProt ID         | <a href="#">MMAD_HUMAN</a>                                                                                                                                                                                                                                              |
| Immunogen Region   | 39-296aa                                                                                                                                                                                                                                                                |
| Immunogen Sequence | MGSSHHHHHH SSGLVPRGSH MGSSDESHVA AAPPDICSRT VWPDETMGPF GPQDQRFQLP<br>GNIGFDCHLN GTASQKKSLV HKTLPDVLAEL PLSSERHEFV MAQYVNEFQG NDAPVEQEIN SAETYFESAR<br>VECAIQTCP E LLRKDFESLF PEVANGKLM I LTVTQKTKND MTWVSEEVEI EREVLLLEKFI NGAKEICYAL<br>RAEGYWADFI DPSSGLAFFG PYTNNTLF |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



MMADHC recombinant protein SDS-PAGE gel  
 3µg reducing conditions Coomassie stain