

Human MMP-10 protein (Recombinant) (His-Tag) (STJP018636)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Proteins |
| Applications | SDS-PAGE/Denatured |
| Host / Source | E.coli |

PRODUCT PROPERTIES


| | |
|---------------------|---|
| Concentration | 0.5 mg/mL |
| Formulation | Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol |
| Molecular Weight | 45.4 kDa |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|--------------------|--|
| Gene ID | 4319 |
| Gene Symbol | MMP10 |
| UniProt ID | MMP10_HUMAN |
| Immunogen Region | 99-476aa |
| Immunogen Sequence | MGSSHHHHH S SGLVPRGSH MGSFSSFPGM PKWRKTHLTY RIVNYTPDLP RDAVDSAIEK ALKVWEEVTP LTF SRLYEGE ADIMISFAVK EHGDFYSFDG PGHSLAHAYP PGPGLYGDIH FDDDEKWTEG ASGTNLFVA AHELGHSLGL FHSANTEALM YPLYNSFTEL AQFRLSQDDV NGIQSLYGPP PASTEEPLVP TKSVP SGSEM PAKCDPALS F DAISTLRG |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human MMP-10 protein under reducing conditions Coomassie blue stain