


Human p130Cas protein (Recombinant) (His-Tag) (STJP018637)

GENERAL INFORMATION	
Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli
PRODUCT PROPERTIES	
Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	43.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	9564
Gene Symbol	BCAR1
UniProt ID	BCAR1_HUMAN
Immunogen Region	465-848aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSRLQQGVS ATVAHLLDLA GSAGATGSWR SPSEPQEPLV QDLQAAVAAV QSAVHELLEF ARSAVGNAAH TSDRALHAKL SRQLQKMEDV HQTIVAHGQA LDAGRGGSGA TLEDLRLVA CSRAVPEDAK QLASFLHGNA SLLFRRTKAT APGPEGGGTL HPNPTDKTSS IQSRPLPSPF KFTSQDSPDG QYENSEGGWM EDYDYVHL
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant human p130Cas protein SDS-PAGE gel under reducing conditions stained with Coomassie blue