

Human MMP-3 protein (Recombinant) (His-Tag) (STJP018675)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol
Molecular Weight	45.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4314
Gene Symbol	MMP3
UniProt ID	MMP3_HUMAN
Immunogen Region	100-477aa
Immunogen Sequence	MGSSHHHHHH SSSLVPRGSH MGSFRITFPGI PKWRKTHLTY RIVNYTPDLP KDAVDSAVEK ALKWVEEVTP LTFSRLYEGE ADIMISFAVR EHGDFYPFDG PGNVLAHAYA PGPGINGDAH FDDDEQWTKD TTGTNLFVA AHEIGHSLGL FHSANTEALM YPLYHSLTDL TRFRLSQDDI NGIQSLYGPP PDSPETPLVP TEPVPEPGT PANCDPALSF DAVSTLRG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing 3µg recombinant human MMP-3 protein under reducing conditions with Coomassie blue staining