

Human HLA-DRA protein (Recombinant) (His-Tag) (STJP018686)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	24.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3122
Gene Symbol	HLA-DRA
UniProt ID	DRA_HUMAN
Immunogen Region	26-216aa
Immunogen Sequence	IKEEHVIIQ AEFYLNPDQS GEFMDFDGD EIFHVDMAKK ETVWRLEEFGRFASFEAQGA LANIAVDKAN LEIMTKRSNY TPITNVPPEV TVLTNSPVEL REPVNLICFI DKFTPPVVNV TWLRNGKPVT TGVSETVFLP REDHLFRKFH YLPFLPSTED VYDCRVEHWG LDEPLKHWE FDAPSPLPET TE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human HLA-DRA protein under reducing conditions Coomassie blue stain