

Human DOK4 protein (Recombinant) (His-Tag) (STJP018694)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M uREA, 10% Glycerol
Molecular Weight	39.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55715
Gene Symbol	DOK4
UniProt ID	DOK4_HUMAN
Immunogen Region	1-326aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMATNFSD IVKQGYVKMK SRKLGIIYRRC WLIVFRKSSSK GPQRLEKYPD EKSVCVLRGCP KVTETISNVKC VTRLPKETKR QAVAIIFTDD SARTFTCDSE LEAEEWYKTL SVECLGSRLN DISLGEPELL APGVQCEQTD RFNVFLLPCP NLDVYGECKL QITHEIYLW DIHNPRVKLV SWPLCSLRRY GRDATRFTFE AGRMCDAG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human DOK4 protein 3µg reducing conditions Coomassie blue stain