

Human RNF114 protein (Recombinant) (His-Tag) (STJP018700)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M uREA, 10% Glycerol
Molecular Weight	28.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55905
Gene Symbol	RNF114
UniProt ID	RN114_HUMAN
Immunogen Region	1-228aa
Immunogen Sequence	MAAQQRD CGGAAQLAGP AAEADPLGRF TCPVCLEVYE KPVQVPCGHV FCSACLQECL KPKKPVCGVC RSALAPGVRA VELERQIEST ETSCHGCRKN FFLSKIRSHV ATCSKYQNYI MEGVKATIKD ASLQPRNVPN RYTFPCPYCP EKNFDQEGLV EHCKLFHSTD TKSVCPCICA SMPWGDPNYR SANFREHIQR RHRFSY

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RNF114 protein under reducing conditions with Coomassie blue stain