

Human APRIL/TNFSF13 protein (Recombinant) (T7-Tag) (STJP018701)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.4M Urea
Molecular Weight	17.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8741
Gene Symbol	TNFSF13
UniProt ID	TNF13_HUMAN
Immunogen Region	105-247aa
Immunogen Sequence	AVLT QKQKKQHSVL HLPINATSK DDSDVTEVMW QPALRRGRGL QAQGYGVRIQ DAGVYLLYSQ VLFQDVTFTM GQVVSREGQG RQETLFR CIR SMPSHPD RAY NSCYSAGVFH LHQGDILSVI IPRARAKLNL SPHGTFLGL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human APRIL/TNFSF13 protein under reducing conditions stained with Coomassie blue