

Human Casein kinase 1 alpha/CSNK1A1 protein (Recombinant) (His-Tag) (STJP018709)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	41 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1452
Gene Symbol	CSNK1A1
UniProt ID	KC1A_HUMAN
Immunogen Region	1-337aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MASSSGSKAE FIVGGKYKLV RKIGSGSFGD IYLAINITNG EEVAVKLESQ KARHPQLLYE SKLYKILQGG VGIPHIRWYG QEKDYNVLVM DLLGPSLEDL FNFCRRFTM KTVLMLADQM ISRIEYVHTK NFIHRDIKPD NFLMGIGRHC NKLFLIDFGL AKKYRDNRT R QHIPYREDKN LTGTARYASI NAHLGIEQSR RDDMESLG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human CSNK1A1 protein under reducing conditions stained with Coomassie blue