

Human LIPG protein (Recombinant) (His-Tag) (STJP018723)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.4M urea
Molecular Weight	38.0 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9388
Gene Symbol	LIPG
UniProt ID	LIPG_HUMAN
Immunogen Region	21-340aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSSPVFPGP EGRLEDKLHK PKATQTEVKP SVRFNLRISK DPEHEGCVLS VGHSQPLEDC SFNMTAKTFF IHHGWTMSGI FENWLHKLVLS ALHTREKDN VVVVDWLPLA HQLYTDAVNN TRVVGHSIAR MLDWLQEKDD FSLGNVHLIG YSLGAHVAGY AGNFVKGTVG RITGLDPAGP MFEGADIIHKR LSPDDADFVD VLHTYTRS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human LIPG protein 3µg under reducing conditions Coomassie blue stain