

Human HOXA9 protein (Recombinant) (His-Tag) (STJP018763)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	32.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3205
Gene Symbol	HOXA9
UniProt ID	HXA9_HUMAN
Immunogen Region	1-272aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MGSMATTGAL GNYVDSFLL GADAADELSV GRYAPGTLGQ PPRQAATLAE HPDFSPCSFQ SKATVFGASW NPVHAAGANA VPAAVYHHHH HHPYVHPQAP VAAAAPDGRY MRSWLEPTPG ALSFAGLPSS RPYGIKPEPL SARRGDCPTL DHTLTLTDY ACGSPPVORE KQPSEGAFFSE NNAENESGGD KPPIDPNNPA ANWLHARS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



15% SDS-PAGE (3µg)

SDS-PAGE gel showing recombinant human HOXA9 protein 3µg under reducing conditions
 Coomassie blue stain