

Human RPL18A protein (Recombinant) (His-Tag) (STJP018766)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	23.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6142
Gene Symbol	RPL18A
UniProt ID	RL18A_HUMAN
Immunogen Region	1-176aa
Immunogen Sequence	MKASGTL REYKVVGRCL PTPKCHTPPL YRMRIAPNH VVAKSRFWYF VSQLKKMKKS SGEIVYCGQV FEKSPLRVKN FGIWLRYSR SGTHNMYREY RDLTTAGAVT QCYRDMGARH RARAHSIQIM KVEEIAASKC RRPVAVKQFHD SKIKFPLPHR VLRRQHKPRF TTKRPNTFF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant human RPL18A protein SDS-PAGE gel 3µg reducing conditions Coomassie blue stain