

Human CD-M6PR protein (Recombinant) (His-Tag) (STJP018771)

GENERAL INFORMATION

| | |
|---------------|----------|
| Product Type | Proteins |
| Applications | SDS-PAGE |
| Host / Source | E.coli |

PRODUCT PROPERTIES


| | |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Concentration | 0.5 mg/mL |
| Formulation | Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT |
| Molecular Weight | 20.3 kDa |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID | 4074 |
| Gene Symbol | M6PR |
| UniProt ID | MPRD_HUMAN |
| Immunogen Region | 27-185aa |
| Immunogen Sequence | TEEKTCD LVGEKKGKESKELALVKRLK PLFNKSFEST VQGGSPTYIY IFRVCREAGN HTSGAGLVQI NKSNGKETVV GRLNETHIFN GSNWIMLIYK GGDEYDNHCG KEQRRRAVVMV SCNRHTLADN FNPVSEERGK VQDCFYLFEM DSSLACSPEI SH |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human CD-M6PR protein 3ug under reducing conditions
Coomassie blue stain