

Human DNMT3L protein (Recombinant) (His-Tag) (STJP018792)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	46.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29947
Gene Symbol	DNMT3L
UniProt ID	DNM3L_HUMAN
Immunogen Region	1-386aa
Immunogen Sequence	MAAIP ALDPEAEPSM DVILVGSSEL SSSVSPGTGR DLIAYEVKAN QRNIEDICIC CGSLQVHTQH PLFEGGICAP CKDKFLDALF LYDDDGYSY CSICCSGETL LICGNPDCTR CYCFECVDSL VGPQTSGKVH AMSNWVCYLC LPSSRSGLLQ RRRKWRSQK AFYDRESENP LEMFETVPVW RRQPVRVLSL FEDIKK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human DNMT3L protein under reducing conditions with Coomassie blue stain