

Human EXOSC3 protein (Recombinant) (His-Tag) (STJP018799)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 0.15M NaCl, 30% Glycerol, 1mM DTT
Molecular Weight	32.0 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51010
Gene Symbol	EXOSC3
UniProt ID	EXOS3_HUMAN
Immunogen Region	1-275aa
Immunogen Sequence	MAEPASV AAESLAGSRA RAARTVLGQV VLPGEELLLP EQEDAEGPGG AVERPLSLNA RACSRVRVVC GPGLRRCGDR LLVTKCGRLR HKEPGSGSGG GYVWVDSQQK RYVPVKGDHV IGIVTAKSGD IFKVDVGGSE PASLSYLSFE GATKRNRPNV QVGDLIYGQF VVANKDMEPE MVCIDSCGRA NGMGVIGQDG LLFKVT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human EXOSC3 protein under reducing conditions Coomassie blue stain