

Human EXOSC1 protein (Recombinant) (His-Tag) (STJP018819)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	23.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51013
Gene Symbol	EXOSC1
UniProt ID	EXOS1_HUMAN
Immunogen Region	1-195aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMAPPVRY CIPGERLCNL EEGSPGSGTY TRHG YIFSSL AGCLMKSSSE GALPVVSVVR ETESQLLPDV GAIVTCKVSS INSRFAKVHI LYVGSMP LKN SFRGTIRKED VRATEKDKVE IYKSF R PGDI VLAKVISLGD AQS NYLLTTA ENELGVVVAH SESGIQMVPI SWCEMQCPKT HTKEFRKVAR VQPEFLQT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human EXOSC1 protein under reducing conditions