

HCV NS5a protein (Recombinant) (His-Tag) (STJP018835)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 40% Glycerol, 1mM DTT
Molecular Weight	13.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	2210-2306aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MSLKATCTAN HDSPDAELIE ANLLWRQEMG GNITRVESENKVVILDSFDP LVAEEDEREV SVP AEILRKS RRFARALPVW ARPDYNPLV ETWKKPDY

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing recombinant HCV NS5a protein 3µg under reducing conditions Coomassie blue stain